DC A-Silicone 5:1

Putty Soft

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DC A-Silicone 5:1

Putty Soft

Elastomeric precision impression material on polyvinyl siloxane basis, addition curing.

DC Dental Central
Grosshandels-GmbH
Carl-Zeiss-Str. 2
D- 22946 Trittau

Made in Germany

GB Instructions for use

Thank you for your confidence in us. With the DC A-Silicone putty soft 5:1 impression material you have made the right decision. Read these instructions for use carefully prior to use. The product is intended exclusively for dental use. Responsible and careful working with the medicinal product is assumed.

Product description:

DC A-Silicone putty soft 5:1 is an impression material that meets the highest demands. The DC A-Silicone putty soft 5:1 material has activated thermoformulation which offers the user a variable processing time and thus the greatest processing comfort. It is patientfriendly as it remains in the mouth for only the same short time regardless of when the impression tray is placed in the mouth. The system is rounded off by the highly contrasting colour system and the pleasant mild mint flavour.

Indication:

Crowns and bridges, partial crowns, inlays, onlays.

Impression technique:

DC A-Silicone putty soft 5:1: Preliminary impression material for the putty wash technique (2-step), situation impressions.
DC A-Silicone putty soft 5:1: Basic material for the simultaneous putty wash technique (1-step), situation impressions.

fig. 1

Processing:

The putty soft 5:1 full cartridges are suitable for use in the following mixers: Pentamix®, Pentamix® 2, Plug & Press® Dispenser und MixStar®. Place the 5:1 full cartridge in the device (please follow the manufacturer's instructions on using the device). Before mounting the mixing tip always start the device and let material flow out until material is flowing evenly from both orifices. Now attach the mixing tip to the orifice and check that it is correctly in place. Then push the fixing ring over the mixing tip as far as it will go and turn it to the right until it locks. Start the mixer again and dispense the desired quantity of material into an impression tray. Hold the tray obliquely and press it on gently (see fig. 1.); while dispensing it, leave the mixing tip in the material. The used mixing tip stays on the cartridge and acts as a cap. Prior to the next use, loosen and remove the fixing ring and remove the used mixing tip, check the orifices for the plug that forms in very rare cases, clean if necessary and proceed as usual.

Impression procedure:

When using a suitable tray adhesive, please follow the manufacturer's instructions for use. Place the tray loaded with the DC putty soft into the mouth, press it on briefly and then hold it in situ until it has set fully. To check the degree of hardening always check the material in the mouth. When using the putty wash technique (two-step) the first impression must always be cleaned and dried in order to ensure curing with the putty wash material.

Desinfection:

The impressions can be disinfected with 2 % glutaraldehyde solution. After removing from the mouth, rinse the impression under running water for 15 seconds. Then disinfect the impression by dipping it completely into disinfectant. Always follow the disinfectant manufacturer instructions.

We recommend the product DC Abformdesinfektion.

Electroplating:

The impressions can be electroplated with copper or silver in the usual baths.

Model fabrication:

The impression can be poured after only 60 minutes. There are no other time restrictions. All class 3 and 4 gypsums that comply with ISO 6873 can be used. The impressions can be cast several times.

Working and safety instructions:

The DC A-Silicone putty soft 5:1 must be used only in accordance with the instruction for use. Any other use, which is not in agreement with the instruction for use, is the sole responsibility of the user.

- The original cap of the 5:1 full cartridges must always be kept. If not used for a prolonged period (approx. 1 week), the cartridge must be closed with this cap to prevent possible formation of a plug. Ensure that the stopper is absolutely clean and that no preliminary curing occurs because the original cap is put on incorrectly.

- Solutions used prior to taking the impression such as oral rinses can interfere with the setting reaction and must be removed thoroughly. Retraction sutures containing iron also have a negative effect on setting.

- Contact between the material and certain types of gloves, e.g. latex gloves, should be avoided as the catalyst can be damaged. We recommend vinyl or PE-based gloves.

- Do not leave any impression material residues in the oral cavity.

- Contact with clothing should be avoided as cured silicones are chemically stable and produce marks that cannot be removed.

- Do not use the DC A-Silicone putty 5:1 material with polyether, polysulfides or condensation-curing silicones.

- Storage: store between 5° C and 27° C.

- Shelf life: 2 years from date of manufacture; see date printed on label. Do not use after the expiry date.

- Processing temperature: room temperature between 18° C and 25° C.

Warnings:

- **Eye contact:** avoid eye contact with the impression material, if necessary, irrigate the eyes immediately with plenty of clear water and consult an ophthalmologist promptly.

- **Swallowing:** after swallowing, rinse out the mouth and drink plenty of water. Always contact a doctor.

- **Allergic reactions:** allergic reactions can occur in sensitive person. In case of doubt, consult a dermatologist or allergologist. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately.

Please note the warning in the safety data sheet also.

Errors and changes reserved.

Technical data and product properties

Technical data	DC A-Silicone putty soft 5 :1
Consistency / EN ISO 4823	Kneadable / Typ 0
Colour	Blue
Hardness Shore A	60
Variable working time including mixing time*	up to 2 min.
Time in mouth	2 min.
Total setting time	up to 4 min.
Deformation under pressure	2,5 %
Recovery from deformation	99,5 %
Linear dimensional change	0,2 %
Flavour	mint

* at 23° C / 73° F. The time specifications in the upper table relate to a relative humidity of 50 +/- 10%. In general, higher temperatures accelerate while lower temperatures delay the setting.

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