D Gebrauchsinformation

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GB Instructions for use

DC Amalgam is a dispersed phase, non gamma-2, amalgam alloy with a high silver concentration. A mixture of lathe cut particles and microfine eutectic spherical particles are packed in self-activating capsules to be used for dental fillings. The handling characteristics are comparable to a conventional amalgam.

Alloy Composition

Silver (Ag) : 69.3%
Tin (Sn) : 19.4%
Copper (Cu) : 10.9%
Zinc (Zn) : 0.4%

DC Amalgam is in compliance of ISO 24234 standards.

Indication:

Stress-bearing fillings in posterior teeth (class I and II).

Contraindications:

- Use of amalgam is not recommended when a physician has deemed it unsuitable for the patient.
- Hypersensitivity (allergic reaction) to amalgam components.

In such cases the use of suitable alternative material must be considered.

Warnings:

Please follow rules and regulations in storing and disposing of dental amalgam and its waste.

- After placement or removal of dental amalgam restorations increased mercury concentration in blood and urine has been observed. According to available scientific knowledge this increase has not been associated with any adverse health effects.

- If placed in close contact with other metal restorations, galvanic effects may occur. In most cases they will be of short lasting duration. If the effect persists, the user should consider replacement of the dental amalgam filling with another material.

- There are no proven adverse effects on the fetus associated with the placement or presence of dental amalgam fillings in the mother. It is sensible however, where clinically feasible, to minimise health interventions during pregnancy and avoid any unnecessary chemical exposure to the fetus. This precaution should be observed with the use of all dental materials.

- In individual cases, local mucosal reactions (lichenoid) have been observed. Such local reactions may be of an irritative (mechanical, chemical, electrochemical) or allergic nature. In case of allergy to components of dental amalgam the use of suitable alternative materials must be considered.

- Unnecessary exposure to mercury vapour or dental amalgam particles during handling, placement or removal of dental amalgam should be avoided. Placement and removal of dental amalgam fillings should be performed with appropriate water spray and vacuum suction.

- Instructions and regulations for storage and disposal of dental amalgam waste must be observed.

For dental use only!

| MIXING TIME IN HIGH ENERGY MIXER | | | | | | | |
|----------------------------------|---------------------------------|---------------------------------|-------------|-------------------|--------|---------|---------|
| Caps, 50 St. | Alloy | Mercury | Mixing time | Revolutions r/min | No. 1 | No. 2 | No. 3 |
| Nr.1 / Yellow REF 716331 | 400 mg Net weight: 20,0 g | 430 mg Net weight: 21,5 g | | 3500 r/min | 9 sec. | 11 sec. | 11 sec. |
| Nr.2 / Orange REF 716332 | 650 mg Net weight: 32,5 g | 700 mg Net weight: 35,0 g | | 4000 r/min | 6 sec. | 8 sec. | 8 sec. |
| Nr.3 / Brown REF 716333 | 800 mg Net weight: 40,0 g | 860 mg Net weight: 43,0 g | | 4500 r/min | 5 sec. | 7 sec. | 7 sec. |


These directions work for the Silamat S5, Ultramat 2, Espe Capmix, and other mixers.

INSTRUCTIONS FOR USE:

Mixing

Please note that the working time and the curing time can vary, depending on the mixer that is used. The mixed amalgam must have a coherent mass and build a shiny surface.

As soon as the mixer is turned on, the bond between the alloy and the mercury occurs automatically.

After the capsule is mixed, open the capsule with the included stick.

Insertion and condensation:

After mixing, the amalgam should immediately be inserted into the cavity using an amalgam carrier. The filling should be built up in layers under adequate condensation pressure using a plugger appropriate for the size of the cavity. Do not use an ultrasonic condenser. Carefully remove excess mercury.

Polishing and finishing:

Approx. 8 - 10 minutes after mixing, the filling can be finished. Polishing should be done in a subsequent session.

THIS ALLOY CONTAINS ZINC; THE AMALGAM MADE MAY SHOW EXCESSIVE EXPANSION IF MOISTURE IS INTRODUCED DURING MIXING OR COMPACTING.

MERCURY TOXIC

- Toxic vapors - avoid inhalation.
- Danger of cumulative effects.
- Keep locked up and out of the reach of children.
- Keep container tightly closed.
- In case of accident or if you feel ill, seek medical advice immediately (show the label when possible).

Recommendations for safety precautions:

Dental clinical personnel, when using silver amalgam, should exercise proper procedures to avoid personal contact with mercury. Precautions should be taken to avoid exposure to mercury vapour in the dental environment. Waste amalgam should be collected for disposal in a manner, which will protect both: those who handle the waste and the environment.

- Treat mercury and amalgam with care.
- Handling of mercury and mixing of amalgam should be done on a separate section of the working table (not close to a source of heat), with raised edges on all sides. This will prevent any spilled mercury from spreading and will facilitate its recovery. The surgery floor should be smooth, without cracks and the surgery and storage area well ventilated.
- Ideally a high-energy mixer should be equipped with a protective cap to cover the capsule while mixing.
- Small droplets of spilled mercury can be collected with the aid of a „mercury-collecting forceps“ or a small amount of freshly mixed amalgam that will easily absorb liquid mercury. NEVER use a vacuum cleaner!
- Placement and removal of amalgam fillings must be done with effective vacuum suction. The use of eye protection and mouth mask is advisable.
- The wastewater drain of the spittoon should be equipped with an amalgam separator. Local regulations on type and installation of such a separator and handling/recycling of the amalgam waste should be taken into consideration.
- Excess mixed amalgam should be kept in a container, designed for the purpose, or at least in a closed container filled with photographic fixer.
- Storage at temperatures no higher than 25° C is recommended.
- Expiry date see bottom.

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control and, therefore, entirely your own responsibility. Should, in spite of this liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

| TEMPS DE MÉLANGE DANS UN MÉLANGEUR À HAUTE ÉNERGIE | | | | | | | |
|--|-------------------------------|-------------------------------|-------------------|-----------------|--------|---------|---------|
| Capsules, 50 St. | Alliage | Mercur | Temps de mixation | Rotations r/min | Nr. 1 | Nr. 2 | Nr. 3 |
| Nr.1 / Jaune REF: 716331 | 400 mg Net contenu: 20,0 g | 430 mg Net contenu: 21,5 g | | 3500 r/min | 9 Sec. | 11 Sec. | 11 Sec. |
| Nr.2 / Orange REF: 716332 | 650 mg Net contenu: 32,5 g | 700 mg Net contenu: 35,0 g | | 4000 r/min | 6 Sec. | 8 Sec. | 8 Sec. |
| Nr.3 / Brun REF: 716333 | 800 mg Net contenu: 40,0 g | 860 mg Net contenu: 43,0 g | | 4500 r/min | 5 Sec. | 7 Sec. | 7 Sec. |

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