## D Gebrauchsinformation

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#### Characteristics

Hardness after casting	ca. 280 HV 10
Specific gravity	8.4 g/cm <sup>3</sup>
Solidus/Liquidus	1,405-1,420° C
Casting temperature	1,470° C
WAK (20° C - 600° C)	14.8 µm/mK
Delivery form	cylinder ø 9.5 x 11 mm
Packaging unit	2.2 lb (1,000 g)

DC bonding alloy NE is a nickel- and beryllium-free cobalt/chrome alloy without carbon.

DIN EN ISO 22674:2007, Type 5

#### Modelling and Sprue System:

To guarantee a safe flow of the modelling, the thickness of the cap should not fall below 0.3 - 0.4 mm. The sprue system of the objects happens as usual with sprue reservoir or crossbar (crossbar should be divided for large bridges). Massive bridge parts should be provided with air outlets.

#### Investing and Preheating:

Any commercial phosphate bonded investments, which are suitable for a preheating temperature of 900° C, are qualified for the investing (please follow the handling instructions of the investment-producer). Large casting objects should only be heated linear.

#### Casting

The optimum preheating temperature of the muffel is 900° C. A ceramic crucible must be used for the melting of the alloy (no graphite crucible). Casting steam should be exhausted.

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**Melting of the alloy with open flame** (acetylene/oxygen or propane/oxygen) with neutral flame without melting powder. When the casting cubes are completely melted, (melting moves under flame pressure), the casting procedure should be started immediately. The oxide layer may not crack, because otherwise alloy components might burn, and die change the characteristics of the alloy!

**When melting the alloy in the high frequency centrifuge**, the casting procedure will be started immediately after all cubes have collapsed and as soon as the shadow over the melting is disappeared. Even here the oxide layer may not crack!

**Cool muffel slowly in the air.**

You should always use new casting cubes, as after repeated casting of the alloy important elements will evaporate and an optimum bonding between metal and ceramic is not longer guaranteed. The elaboration of the modelling happens with suitable milling cutters and Al-Oxyd stones. The thickness of the cap should not go below 0,2 - 0,3 mm after finalization

#### Firing of the ceramic

An oxide firing is not necessary. If an oxide firing is desired, then under atmosphere (wash brand is not applicable for bonder brand). Then sandblast with aluminum oxide 110µm. Never put non-precious dental alloys in acid-stain bath.

#### Soldering and Welding

Commercially available solder on cobalt base is suitable for soldering. Never use gold- or palladium solder for the soldering of NE objects among themselves.

#### Cleaning

In ultrasonic bath or with steam blasting.

**Guarantee:** All recommendations regarding the handling are based on the manufacturer's experiences and tests. Therefore they may only be taken as guide values. The products are subject to a continuous advancement. Construction and composition are subject to alteration.

#### Risks phrases:

R42/43: (Kobalt) May cause sensitization by inhalation and skin contact.  
R53: (Kobalt) May cause long-term adverse effects in the aquatic environment.

#### Safety phrases:

S2: Keep out of the reach of children.  
S22: Do not breathe dust. (Kobalt)  
S24: Avoid contact with skin. (Kobalt)  
S37: Wear suitable protective clothing and gloves. (Kobalt)

#### Disposal:

Follow the disposal considerations in the material safety data sheet

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