DC Bite & DC Bite extra hart

Stand 06 / 2012

DC DentalCentral GmbH • Owidenfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

D Gebrauchsinformation

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GB Instructions for use

Polyvinyl siloxane bite registration materials – DIN 13903.

DC Bite und DC Bite extra hart are syringeable addition-curing elastomeric bite registration materials intended for making accurate occlusal records.

Indications

DC Bite and DC Bite extra hart can be used for making accurate occlusal registrations.

DC Bite extra hart can be used for applications, which require a hard bite registration silicone.

Contraindications

Do not use DC Bite und DC Bite extra hart for conventional impressions as their high end hardness values severely inhibit withdrawal from the mouth and may result in inadvertent extraction of teeth.

Cautions

- Do not use after expiration date.
- Do not leave residual bite registration material in the sulcus or oral cavity.
- Do not swallow bite registration material! If swallowed: consult a medical doctor if problems arise or persist.
- Avoid contact with the eyes. If there is accidental contact with the eyes: rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.
- The product does not normally cause allergic reactions; however, for sensitive persons, an allergic reaction to the product cannot be ruled out.
- Use of products containing active sulfur, aluminum chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, latex gloves, etc.) in conjunction with this product will interfere with the setting reaction of the polyvinyl siloxane material. Use of these materials requires the area to be rinsed thoroughly to remove all residues.
- If gloves are worn when handling the material, they should be tested with a sample of the bite registration material beforehand to ensure compatibility. Avoid contact with latex gloves.
- Do not allow the material to enter the sewerage or water system to avoid environmental hazard.
- Avoid contact with clothing, as the material cannot be removed by dry cleaning.

Note

The total working time and intraoral setting time should be strictly adhered to in order to ensure optimal bite registrations. Should the temperature of the material deviate noticeably from 23° C (74° F), the total working time and setting times will be affected.

Instructions for Use

- Insert the cartridge into the DC safety dispensing gun S50 1:1/2:1. Ensure that the notches on the cartridge base are pointing downwards. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.
- Remove the cartridge cap. The cap can be replaced after initial use.
- Advance plunger into cartridge and dispense a small amount of bite registration material until equal amounts are extruded at the same rate.
- Install a mixing tip on the cartridge, turn cap 90 degrees in a clockwise direction to lock in place.
- Inject material directly onto the patient's occlusal surfaces or occlusal rims. Leave used mixing tip on the cartridge after use or replace the cartridge cap.
- Before using the cartridge again, remove cartridge cap or remove and discard the used mixing tip. Check the cartridge openings to ensure no polymerized material is present. Should this be the case, unblock the openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.
- Install a new mixing tip and continue as described above.

The DC safety dispensing gun S50 can be sterilized in the autoclave.

Storage Conditions:

In a dry place, at 18 – 25° C (64 – 77° F). Protect from light.

Disinfection

Bite registrations can be disinfected using, e. g. 2 % acidic glutaraldehyde solution. See manufacturer instructions.

The information provided for DC products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Product	Total working time at 23° C / 74° F (minimum)	Intraoral setting time at 35° C / 95° F (minimum)	Total setting time* (minimum)	Hardness (minimum)	Linear dimensional change (maximum)	Mixing technique	Mixing element
DC Bite	30 sec.	1 min. 30 sec.	2 min.	Shore-A 90 / HD 58	-0,20 %	DC safety dispensing gun S50	DC mixing tip, ø 6,5 mm, REF 701243
DC Bite extra hard	30 sec	1 min. 30 sec.	2 min.	Shore-D 43 / HD 76	-0,10 %	DC safety dispensing gun S50	DC mixing tip, ø 6,5 mm, REF 701243

*Total setting time (removal time from mouth) from start of mixing.

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