DC Knetmasse & DC C-Silicone

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Straße 2 • D-22946 Trittau

D

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DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Straße 2 • D-22946 Trittau

GB

Silicon based impression material, condensation curing.

The DC system compliments all normal applications or can be combined individually.

Product	ISO 4823	Colour	DOSAGE			Mixing-time	Total working time	Time in mouth	Setting time	Hardness Shore A	Strain in compression	Recovery from deformation	Linear dimensional change
			Basis	DC Hardenerpaste	DC Hardenerliquid								
DC Putty <i>soft</i>	type 0	white	1 measuring spoon 12 ml = (22,4 g)	2,5 graduation = 2,5 cm 0,53 ml = (0,5 g)	10 drops 0,24 ml = (0,24 g)	30 sec.	1 min. 45 sec.	4 min.	5 min. 45 sec.	60	2 %	98,8 %	ca. 0,2 %
DC Putty <i>hard</i>	type 0	white	1 measuring spoon 12 ml = (21,7 g)	2,5 graduation = 2,5 cm 0,53 ml = (0,5 g)	10 drops 0,24 ml = (0,24 g)	30 sec.	1 min. 45 sec.	4 min.	5 min. 45 sec.	70	2 %	99,4 %	ca. 0,3 %
DC Silicone <i>green</i>	type 1	green	10 graduation 6,1 ml = (8,5 g)	5 graduation = 5 cm 1,0 ml = (1,0 g)	10 drops 0,24 ml = (0,24 g)	30 sec.	1 min. 45 sec.	4 min.	5 min. 45 sec.	50	5 %	99,5 %	ca. 1 %
DC Silicone <i>blue</i>	type 2	bue	10 graduation 5,8 ml = (7,5 g)	5 graduation = 5 cm 1,0 ml = (1,0 g)	10 drops 0,24 ml = (0,24 g)	30 sec.	1 min. 45 sec.	4 min.	5 min. 45 sec.	40	6 %	99,7 %	ca. 1 %

1 graduation = 1,0 cm (related to the mixing pad)

Results from 23° C, 50 ± 10 % rel. humidity, higher temperatures shorten, lower temperatures prolong the setting times.

Typ 0 = putty consistency – initial impression for the correction material or base material for the double mixing technique

Typ 1 = heavy-bodied consistency – mono phase impression

Typ 2 = medium-bodied consistency – correction material for the double impression and double mixing technique

Preparation and processing:

The contrasting colours of the base and catalyst masses ensure good control during mixing. Putty material and hardener needs to be carried out manually and the correction material needs to be carried out manually with a spatula according to the dosage.

Carefully close tubes and containers after use.

Disinfection: Impression may be disinfected with 2 % glutaraldehyd-solution. We recommend DC impression disinfection.

Electroplating: The silicone material is compatible with the customary dental acid or, respectively, alkaline copper or silver electroplating baths.

Modell preparation: Pouring of the impression should be carried out no sooner than 60 minutes after removal of the impression material from the mouth; and pouring up to 7 days possible, e.g. with type 3 or 4 plasters that meet the requirements of ISO 6873.

Storage: Do not store above 27° C.

Shelf life: see information on package

Handling and safety measures:

Avoid contact with eyes, if contact should occur, rinse eyes immediately and thoroughly in clean running water and go directly to an eye specialist. Avoid getting traces of the impression mass on clothing, as its chemically bonds to clothing fibers. Close paste hardener vessels carefully after use.

Warnings for DC Hardenerpaste / DC Hardenerliquid:

Inhalation: bring the person into the fresh air and treat with oxygen if necessary.

Eye contact: Avoid eye contact, if necessary, irrigate the eyes immediately with plenty of clear water and consult an ophthalmologist promptly.

Skin contact: Avoid prolonged skin contact. After prolonged contact wash the skin with soap and water.

We recommend wearin suitable protecting clothing and gloves.

Swallowing: after swallowing, rinse out the mouth and drink plenty of water. Always contact a doctor.

Allergic reactions: allergic reactions can occur in sensitive person. In case of doubt, consult a dermatologist or allergologist. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately.

The preparation, usage and application of the products supplied lie exclusively in your area of responsibility since we have no opportunity of exercising control over these matters. The advice given on technical application matters verbally, in writing or through trials, then this advice is to be understood as nonbinding information; the same applies in respect of the possible protective rights of third parties. Such advice does not release you from the obligation of testing the products supplied for their suitability for the intended purpose and processes. Our liability is an every case of damage limited to the value of the materials supplied and the goods used. The proper quality of our products is guaranteed.

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