DC DON hot

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Str. 2 • D-22946 Trittau • www.dental-central.de

D Gebrauchsinformation

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GB Instruction for use

Thank you for choosing DC Don hot, the high-tech hot-curing acrylic.

DC Don hot allows uncomplicated handling and is, beside the packing- or pressing technique, also perfectly suitable for the injection technique. To make sure you always achieve reliable results with this product we would like to ask you to read these instructions thoroughly.

Indications:

- Total upper and lower dentures by using the pressing- / or packing technique
- Total upper and lower dentures by using the injection technique

Short Information:

Mixing ratio powder/liquid:	10 : 4
Casting phase approx.	10 min*
Processing- & pressing phase approx.	45 min*
Polymerization at 95° C	20 min*

* The provided processing times refer to a material- and room temperature of 23° C [73,5° F]. Higher temperatures extend, lower temperatures shorten these processing times.

About this product:

DC Don hot pink belongs to the newest generation of short-term hot-curing dental acrylics. No complex equipment is required. The comfortable plastic consistency assures a very satisfying performance during the pressing process and guarantees an accurate fitting denture.

Material:

Denture base resin, short-term hot-curing acrylic based on methyl-methacrylate in the form of powder and liquid. Classification according to EN ISO 20795-1 Type 1 Group 1 and according to medical device directive 93/42/EEC annex IX Class IIa for removable and for fixed dentures.

Product features:

- very good flow ability
- high mechanical properties
- physiological shades and 100 % color stable
- high quality results and an accurate fitting denture
- easy handling
- mucous membrane tolerance due to a low amount of residual monomer
- compatible to all known individual shades
- evaluated and certified biocompatibility

DC Don hot does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded) and/or supplied with retentions.

Mixing Ratio 10 : 4

Suggested mixing ratio: 10 g of powder with 4 ml of liquid. An individual dosage is possible. For an extended processing phase it is possible to mix 10 g of powder with 5 ml of liquid. By doing that the dough can be handled for 2 - 3 hours without affecting the high mechanical properties.

Preparation:

The wax modellation is embedded into the stable plaster. The flask is boiled out thoroughly to remove the wax before applying alginate based insulation on the surface of the plaster. Important: The insulation needs to be well-dried to avoid white marks or spots on the finished acrylic.

Processing:

Place required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as required to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a thickening phase of approximately 5 to 10 min. the material can be handled within another 45 min.

Pressing:

Place the material with a little excess into the flask. Slowly close the flask and continuously apply pressure. Keep the pressure for at least 10 min. A test pressing is not necessary due to the excellent flow ability of the acrylic.

Polymerization:

This acrylic is especially designed for the short-term hot-curing polymerization. Immediately place the flask in boiling water (approx. 3 l of water per flask) and allow curing for about 20 min.

For curing larger projects place the flask in boiling water, turn heat source off and wait 15 min. Bring water to boil again and allow curing for 10 min.

Of course the material can be processed by using any other polymerization technique as well. Allow flask to cool down as usual.

After completion:

Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste.

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25° C [77° F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapours. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

Number of units:

- 714151 DC Don hot 1.000 g powder
- 714156 DC Don hot 500 ml liquid