DC Dublier Silicone & DC Dublier Silicone Quick

DC DentalCentral GmbH • Owidenfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

D Gebrauchsinformation

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GB Instructions for use

DC Dublier Silicone & DC Dublier Silicone Quick is an addition-curing 1:1 duplicating silicone with final hardness of Shore A 22 and low viscous with best mechanical characteristics. It offers an extremely high accuracy in dimension and design of the duplicating form and the pink tint supports an optimal surface characterisation.

DC Dublier Silicone & DC Dublier Silicone Quick is very flexible in the final consistency despite high ultimate tensile strength and therefore easy to deflask particularly in case of undercuts. It is from extraordinary productivity since it contains no additional fillers and requires no wetting agents. Due to a platinum catalyst in component A the hardening without any shrinkage is ensured.

DC Dublier Silicone & DC Dublier Silicone Quick is also suitable for the casting without flasks and causes no limitation to the expansion value of the investment.

Working instructions

DC Dublier Silicone & DC Dublier Silicone Quick component A and B must be mixed at room temperature in proportion 1:1 approx. 1 min. until an uniform colour is reached. Pour it from a height of approx. 30 cm constantly onto the dry plaster model in the flask until it's all round covered. If applicable in case of manual mixing removal from bubble formation should be indicated by a vibrator.

After a curing period of approx. 30 minutes in case of DC Dublier Silicone or 10 minutes in case of DC Dublier Silicone Quick at a room temperature of 23° C the model can be replaced and should be left another few minutes to release probably existing reproductions of undercuts.

For further questions don't hesitate to contact us.

Shelf life

DC Dublier Silicone is storable in original locked containers at least 12 months at a temperature between 10° C and 25° C.

uncured	DC Dublier Silicone & DC Dublier Silicone Quick	
	Component A	Component B
Viscosity	2.500 mPa s	2.500 mPa s
Density (20° nach DIN 51 757)	1,075 g/cm ³	1,075 g/cm ³
Colour	white	pink

cured	DC Dublier Silicone	DC Dublier Silicone Quick
	Processing time	6 minutes
Setting time	30 minutes	10 minutes
Shore A (DIN 53 505)	16-22 after 30 minutes	16-22 after 30minutes
	18-22 after 60 minutes	18-22 after 60 minutes
Tear resistance (DIN 53 504 S1)	3,5 N/mm ²	3,5 N/mm ²
Elongation at tear (DIN 53 504 S1)	400 %	400 %
Continuing tear resistance (ASTM D 624 B)	4,5 N/mm ²	4,5 N/mm ²
Colour	pink	pink
Linear change of measurement	< 0,2 %	< 0,2 %

These data represent approximate values and are not intended to produce specifications.

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