## D Gebrauchsinformation

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## GB Instruction for use

Thank you for choosing DC DUR cold, the top acrylic for denture repairs.

DC DUR cold is an economic acrylic for quick and easy repairs of dentures and orthodontic devices. To make sure you always achieve reliable results with this product we would like to ask you to read these instructions thoroughly.

### Indications:

- Repairs of total and partial dentures
- Repairs of orthodontic devices
- Relinings
- Completion of dentures

### Short information:

Mixing ratio powder/liquid:	10 : 6
Swelling phase approx.	30 sec*
Pouring phase approx.	3 min*
Polymerization under pressure at approx. 30° C	10 min*

\* The provided processing times refer to a material- and room temperature of 23° C [73,5° F]. Higher temperatures extend, lower temperatures shorten these processing times.

### About this product:

**No matter if a total denture is fractured or broken in half, or the saddle from a model cast denture is damaged, DC DUR cold can put things right in minutes.**

### Material:

Denture repair resin, self-curing acrylic based on methyl methacrylate in the form of powder and liquid. Classification according to DIN EN ISO 20795-1 Type 2 Group 1 and according to medical device directive 93/42/EEC annex IX class IIa for removable and for fixed dentures.

### Product features:

- very good flow ability
- easy handling
- excellent chemical bonding with any denture base acrylic
- color stable due to UV- and thermal stabilizers
- fast and safe curing

**DC DUR cold does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grounded) and/or supplied with retentions.**

### Mixing ratio 10 : 6

Suggested mixing ratio: 10 g of powder with 6 ml of liquid. For use in casting procedures an individual dosage is possible.

### Preparation:

The edges of the pieces that need to be repaired should be grinded and slightly beveled to ensure a secure bonding. The edges also need to be completely free of wax residue.

### Isolating:

When using a plaster model it is necessary that it is isolated carefully.

Important: The isolation needs to be well-dried to make sure there are no white marks or spots on the finished acrylic.

### Processing:

Place required amount of liquid into the mixing cup, then add the according amount of powder. Make sure a small liquid surface remains on top. Stir thoroughly but keep the material "bubble-free". Apply the still fluid material with a little excess on the surface in need of repair. When thicker layers are needed it is necessary to apply several thinner layers. Each layer has to be completely cured and cooled down before applying the next layer.

Important: The material must not be applied on the surfaces as soon as it reaches its plastic phase. Otherwise a secure bonding of the two materials cannot be guaranteed.

### Polymerization:

The material cures - when processed at room temperature of approx. 23° C - within 10 to 12 minutes. Curing the material in a polymerizing apparatus (2 - 3 bar air pressure, 30 - 40° C water temperature) does largely eliminate the formation of pores or trapped air which occur as a result of exothermal generation of heat. The repaired parts can be grinded or polished immediately after the material is cured.

### After completion:

Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste.

### Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25° C [77° F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

### Number of units:

714141 DC DUR cold 1.000 g powder  
714146 DC DUR cold 1.000 ml liquid