# DC EugenolfreeCem

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Str. 2 • D-22946 Trittau • Germany • www.dental-central.de

Stand 09 / 2010

## D Gebrauchsinformation

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## GB Instructions for use

DC EugenolfreeCem is an eugenol-free easily to process cement. It does not influence the setting of resin-based filling and luting materials.  
It is recommended to use DC EugenolfreeCem in case of allergies to eugenol. To support the preventives of tooth vitality DC EugenolfreeCem contains calcium hydroxide and hydroxylapatite.

**Field of Application**  
DC EugenolfreeCem is used as a temporary cement for attachment and sealing.

**Method of use**  
DC EugenolfreeCem should be used at room temperature. In case of lower temperature (use of a cool glass-plate) the working time will be extended correspondingly. Base and catalyst paste should be mixed in proportion: 2 : 1 v/v.  
For that purpose, put one length of the catalyst paste and accordingly double the length of base on a mixing pad of strong, non-absorbing material or glassplate respectively (1 : 1 w/w). Then mix the two pastes for approximately 30 seconds until the material takes on a homogeneous consistency. Use a plastic spatula.

Apply a thin layer inside the temporary crown or bridge and affix to the prepared tooth. Do not use dry cores, moisten them with some water beforehand. DC EugenolfreeCem remains workable for about 1-2 minutes and sets in the mouth within 3-6 minutes. 3-5 minutes after DC EugenolfreeCem has set remove any excess.

**Supplementary indications**  
To obtain a softer material and to keep DC EugenolfreeCem workable for a longer time, add a little bit more of DC EugenolfreeCem-base or a small trace of Vaseline. A humid milieu (humid cores) additionally facilitates to remove work afterwards. A harder and faster setting material is obtained by adding more catalyst paste (temporary sealing), or zinc oxide. Extensive superficial and multisurfaced cavities as well as temporary longterm sealing work (max. 1 month) should be treated with a harder cement. Residues of DC EugenolfreeCem can be removed from crowns easily after they have been shortly sprayed with a cooling spray.

**Side Effects**  
Undesirable side effects of this medical-product are not to expect if working correspondingly to the instruction of use. In case of side effects – even when doubtful – please inform us with exact description of accompanying circumstances and symptoms.

**Instruction references**  
Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made to the instruction of use) please inform us, if possible, with exact description of accompanying circumstances and symptoms.

**Warning**  
Avoid contact with skin or mucus membranes. If contact does occur, use alcohol-soaked cotton rolls to dissolve the cement. Rinse with plenty of clean water. If eye contact occurs, wash contaminated eye lid with plenty of clean water and seek medical attention.

**Storing**  
Store DC EugenolfreeCem at room temperature. Close tubes tightly after use. Do not exchange the close-ups of the tubes.

**Presentation**  
25 g base in tube,  
25 g catalyst in tube

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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