# DC Feinsthybrid NF

DC DentalCentral GmbH • Owiedenfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

## Gebruuchsinformation

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Entspricht:  
EN ISO 4049

## GB Instruction for use

DC Feinsthybrid NF is a light-curing, radioopaque nanohybrid composite material for inlay therapy. DC Feinsthybrid NF is cured at a wavelength of 400 – 500 nm (blue light).

### Composition:

The fillers consist of barium glass fillers, radio-visible fillers and silicone dioxide with a particle size between 0,04 and 3 µm. The inorganic fillers amount to a portion of 76% (w/w) or 60% (v/v), respectively. DC Feinsthybrid NF comprises methacrylates (22% w/w). Furthermore it comprises additives, initiators, stabilizers and dyes (0,8% w/w).

### Indication:

- DC Feinsthybrid NF is a universal inlay material for direct restoration of cavities of anterior and posterior teeth (also occlusal stress bearing).
- Class III and IV anterior tooth inlays.
- Class V inlays (cervical caries, root caries, wedge-shaped defects).
- Class I and II posterior area inlays.
- Veneering for discoloured anterior teeth.
- Fixation of loose teeth.
- Fissure sealing in the molar and posterior area.
- Repair of composite and ceramics veneers.

### Contraindication:

Existing allergy to components of DC Feinsthybrid NF.

### Interactions:

In rare cases, ingredients may cause sensitisation in particularly susceptible persons. In such cases, DC Feinsthybrid NF should not be used. In order to reduce the risk of tooth pulp reactions, bleeding dentine should be sealed with a suitable material to protect tooth pulp and dentine. Phenolic and eugenol-comprising substances may inhibit polymerisation. They must not get into contact with DC Feinsthybrid NF. Discolourations may occur in combination with cationic mouth washes, plaque indicators and chlorhexidine.

### Application:

1. **Colour selection:** Prior to colour determination, the tooth should be cleaned and polished. The colour is chosen with the tooth moist.
2. **Drying:** Careful drying using rubber dam is mandatory.
3. **Cavity preparation:** Prepare cavities according to the rules of adhesive technology in order to protect, among other things, the tooth substance. Avoid sharp edges and cuts. Do not prepare caries-free cervical defects but only clean them using pumice or other suitable cleaning paste and rubber polishers or cleaning brushes. Then clean the cavity using water spray.

4. **Pulp protection/lining:** For very deep areas near the pulp, introduce a calcium hydroxide preparation, then overlay with pressure-resistant cement (e.g. DC Glass-Ionomercement). The cavity walls should not be covered since they will be needed for bonding with enamel/dentine adhesive (e.g. DC Adhesive).
5. **Dies/interdental wedges:** Use die ribbon if required.
6. **Enamel and dentine conditioning:** Apply 37% etching gel (e.g. DC Ätzel) to enamel and dentine (beginning at the enamel edges). Follow the instructions for the use of the product applied.
7. **Bonding:** As the products by DC Dental Central are designed for compatibility, we recommend to use DC Adhesive. Follow the instructions for the use of the product applied.
8. **Application of DC Feinsthybrid NF:** DC Feinsthybrid NF is introduced at a layer thickness of no more than 2 mm and shaped using a suitable instrument. Each layer is cured with a light lamp for 40 seconds. In doing so, the light conductor is held as close to the material as possible. If metal dies were used, the buccal or lingual/palatal sides, respectively, should be light-cured once more after removal of the die.
9. **Completion of the restoration:** After light-curing, excess material is removed using suitable finishers or fine diamond grinders. Control occlusion and articulation and adapt if necessary. Finally, polish to high gloss using silicone polishers, disks or stripes.

### Further information:

1. If necessary, DC Feinsthybrid NF may be applied directly to cured material. If the inlay has been polished before, it should be roughened first, and bonding should be applied to the surface prior to applying new DC Feinsthybrid NF.
2. DC Feinsthybrid NF should be used at room temperature only (not at temperatures exceeding 28 °C). At lower temperatures, the material is more difficult to squeeze out.

### Caution:

Avoid contact of DC Feinsthybrid NF with skin, eyes and mucous membranes. Unpolymerized material may cause slight irritations and sensitisation to methacrylate. Commercial medical gloves do not provide protection against sensitization to methacrylates. Therefore, contact with gloves should also be avoided.

### Storage:

- Do not use DC Feinsthybrid NF after the stated expiry date.
- Store at 2 – 28 °C.
- After use, immediately seal the syringe.
- Expiry date: see stamp on the syringe/DFT capsule packaging.
- **Keep out of the reach of children!**
- For use in dentistry only!

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