DC Nahtmaterial / DC sutures / DC matériaux de suture

DC DentalCentral GmbH • Owidienfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

D Gebrauchsinformation

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GB Instruction for use

Before using this device read instructions carefully. The product is to be used for the described indications only. The product should be inserted, manipulated and removed only by a qualified, licensed physician or under his direction. The medical techniques described in these instructions do not represent all medically acceptable procedures, nor are they intended to substitute the physician's judgement in treating any specific patient. Do not use a product when the individual package is broken or dirty or if the product itself looks abnormal in any way. Use the device prior to the expiry-date on the package label. Do not sterilize. DC DentalCentral does not warrant for expired or resterilized products.

Technical data:

Technical data and components are indicated on the product label. Further ref. numbers and configurations are obtainable from the product catalogue.

Material:

Needle: ST. Suture: see label

Indications and Applications:

Indications are indicated on the product label.

Contraindications:

Owing to the in vivo decreasing tensile of surgical suture material made of nylon or silk, it should not be used in patients, if a permanent retention of the tensile strength is necessary.

Surgical suture material made of silk is not approved for use in patients who are sensitive towards or allergic to silk.

Surgical suture material made of silk is generally not approved for use in cardiovascular surgery.

Contraindications in terms of the use of surgical suture material made of polypropylene, polyester, Supramid (PA) or Medilene (PVDF) are unknown.

Warnings:

As with foreign bodies, the long-term contact of surgical suture material made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF) with salt solution, e.g. with those in the urinary or gall bladder systems, can lead to the formation of concretions.

Since the risk of wound dehiscence is dependent on the area of the body and suture material being inserted, the user should be familiar with surgical methods and techniques (especially knot techniques), in which surgical suture materials made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF) are used for sealing wounds. In the event of drainage, or sealing of infected or contaminated wounds, generally acknowledged surgical techniques must be followed.

In general, when using surgical suture materials made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF), you should ensure that the suture material is not damaged or squashed by surgical instruments such as forceps, needle-holders or tweezers, since this could result in negative changes in the tensile strength of the suture material and in its safety.

When using surgical suture material made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF), one should ensure to grip the middle of the needle and not to damage either the point or armed end of the needle with surgical instruments like tweezers or needle holders, since otherwise this can cause the needle to bend or break. In general, surgical needles should be handled with meticulous care when in use to avoid needle-prick injuries to patient and user.

Suture material made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF) causes a slight acute inflamed reaction in the tissues. This leads to a general encapsulation of the suture material by fibrous connective tissue.

Surgical suture material made of nylon, silk, polypropylene, polyester, Supramid (PA) or Medilene (PVDF) is non-absorbable.

Advanced hydrolysis of the nylon thread resp. progressive degradation of the perfluorinated silk fibre in vivo can lead to a decrease in the tensile strength of the thread.

Surgical suture material made of polypropylene, polyester, Supramid (PA) or Medilene (PVDF) is non-absorbable. Reduction in the tensile strength in vivo is unknown.

Side-effects

The following side-effects can occur with the insertion of surgical suture materials made of nylon or silk: wound dehiscence, decrease in the tensile strength of the thread in vivo, concretum formation due to long-term contact with salt solution in the body, wound infections, slight acute inflamed reactions in the tissues and temporary local irrigation. In addition, allergic reactions in patients sensitive towards silk can occur.

The following side-effects can occur with the insertion of surgical suture materials made of polypropylene, polyester or Medilene (PVDF): wound dehiscence, concretum formation due to long-term contact with salt solutions in the body, wound infections, slight acute inflamed reactions in the tissues and temporary local irritation.

Sterility:

Suture materials made of nylon, silk, polypropylene, polyester and Medilene (PVDF) are single packed and sterile. Do not use the suture material if the packing is damaged. The sterility for unopened packs is guaranteed until the expiry date.

Storage conditions:

Suture materials made of nylon, silk, polypropylene, polyester and Medilene (PVDF) should be stored at a temperature below 32°. Protect from dampness and direct sunlight. Do not continue using after the expiry date.

Disposal:

After use discard the product according to the regional or national waste regulations for infectious materials.

Exclusion of liability:

We declare on our own responsibility that the product is in compliance with the directives of 93/42/EEC for medical devices. In consideration of the biological differences of patients, the efficiency of the product cannot absolutely be guaranteed. DC DentalCentral cannot guarantee the success of the product, or exclude the occurrence of undesirable effects, as the conditions of diagnosis and indication, as well as of application and use of the product are beyond our control. These instructions have been written in German originally and have been translated into English. In case of discrepancies in interpretation the German version has priority.

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