DC PRESS regular

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Str. 2 • D-22946 Trittau • www.dental-central.de

D Gebrauchsinformation

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GB Instruction for use

Thank you for choosing DC PRESS regular, the universal dental acrylic.

DC PRESS regular is the all-rounder among all denture base acrylics and especially designed for dental technicians who like to cover all indications with only one product. To make sure you always achieve reliable results with this product we would like to ask you to read these instructions thoroughly.

Indications:

- Completion of model cast dentures
- Indirect relinings
- Immediate dentures
- Dilatations and repairs

Short information:

Mixing ratio powder/liquid:	10 : 7
Thickening phase approx.	30 sec*
Pouring phase approx.	3 min*
Plastic phase approx.	3 min*
Pressure polymerization at 45 - 50° C approx.	7 min*

* The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures extend, lower temperatures shorten these processing times.

About this product:

Whether it is the completion of a partial denture, a simple repair, an indirect extraoral relining or a dilatation, DC PRESS regular always shows reliable results while saving precious time.

Material:

Denture base resin, Autopolymer based on methyl methacrylate in the form of powder and liquid. Classification according to DIN EN ISO 20795-1 type 2 group 1 and according to medical device directive 93/42/EEC annex IX class IIa for removable and for fixed dentures.

Product features:

- cadmium free
- very good flow ability
- very high fitting accuracy
- color stable due to a catalyst system without tertiary amines
- mucosal friendly due to a low amount of residual monomer
- fast and safe curing
- certified biocompatibility

DC PRESS regular does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded) and/or supplied with retentions.

Mixing ratio 10 : 7

Suggested mixing ratio: 10 g of powder with 7 ml of liquid. Dosing individually is possible.

Preparation:

The arrangement of the prosthetic teeth in wax need to be fixed with a key made of plaster or silicone. Key, model and teeth are boiled out thoroughly to remove any wax residue.

The plaster model needs to be well watered to assure that no air is pressed from the plaster into the acrylic while curing inside the pressure polymerization unit.

Isolation:

Apply alginate based isolation on the model. Important: The isolation needs to be well dried to make sure there are no white marks or spots on the finished acrylic.

Processing:

Place required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as needed to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a thickening phase of approximately 30 sec. the material is pourable for another 2 to 3 minutes. Within this phase the material should be poured into the pre-casts. As soon as the plastic phase begins the material stops flowing and can be modelled for approx. 3 minutes.

Polymerization:

After a maximum of 8 minutes the polymerization should be carried out in a pressure polymerization unit at a pressure of 2 to 3 bar and a water temperature of 45 to 50° C for 6 to 7 minutes.

After completion:

Until incorporation the denture should be stored in water. It is recommended to store the denture in water for at least 24 hours to achieve an even higher resorption of the residual monomer and thus reduce irritations of taste.

Storage and safety measures:

Material must not be used after expiration date. Storage temperature must not exceed 25° C [77° F]. Avoid direct sunlight. Keep containers closed after use and keep out of reach of children. Monomer is highly flammable. Irritating to respiratory system and skin. Sensitization may occur after skin contact. Keep containers in a well ventilated place. Avoid eye and skin contact. Do not inhale vapours. Keep away from sources of ignition. No smoking. Do not allow to get into drains. Empty containers can be recycled in accordance to your local recycling station.

Number of units:

- 714111 DC PRESS regular 1.000 g powder
- 714116 DC PRESS regular 1.000 ml liquid