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DC DentalCentral GmbH • Owidenfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

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
Item no. 714442 DC Polymerisationsdrucktopf 4 l. without inflating hose

Serial no.: _____

Before start up please read the operating instructions of the equipment attentively.

The DC Polymerisationsdrucktopf is a pressure vessel for polymerisation that is protected with special packing against damage during shipment. If nevertheless a transport damage should be ascertained please inform the forwarder and the sender immediately.

1. Installing of the equipment

Connect the enclosed or existing inflating tube at 3 bar (max. 8 bar) to the compressed air network, if necessary install a pressure-reducing valve between. Do not place the equipment in the immediate vicinity of uncontrolled heat sources, such as open fire. Always comply with Pressure Equipment Directive and general accident prevention instructions when handling the DC Polymerisationsdrucktopf.

2. Handling of the equipment Pour in water

The DC Polymerisationsdrucktopf should be filled with lukewarm water (approx. 40° C / 104° F) to 3-5 cm off the rim.

Put on the lid

Place the lid fully down on the bayonet catches on the vessel body. The marking line "I" on the lid must be aligned with the locking mechanism on the handle. Turn the red switch to the open (AUF) "I" position. Rotate the lid clockwise until both handles are in line. Turn the red switch to the closed (ZU) "0" position. By that the pressure release valve is closed and the lid is mechanically locked. Thus it cannot be opened under pressure.

Build up pressure

Pressurize the vessel with the inflating tube to 2.0-2.2 bar (green marking line) by pressing the plug nipple on the inlet valve (check if the transport protection is twisted off).

Should the pressure exceed the green marking line it can be reduced by briefly opening the red switch. Should the pressure exceed 2.5 bar (red marking line) the safety valve automatically opens and reduces the pressure to 2.0 bar.

Polymerisation

The polymerisation time usually amounts 30 min and the lid should remain closed during this time to ensure a pore-free hardening of the object (Consider the acrylic manufacturers instructions). To insert further objects a briefly opening is possible earliest after 5 min.

Releasing the pressure and opening

Turn the red switch to the open (AUF) "I" position. Wait until pressure has been released completely. Turn the handle of the lid anti-clockwise until unlocked (the marking line on the lid should be aligned with the locking mechanism on the handle from body of the vessel) and then lift off.

3. Maintenance and care of the equipment

Build-up of lime can prevent the safety valve from functioning. The safety valve should be checked twice a month by unscrewing the serrated nut. Additional checking every 6 months:

Fill the DC Polymerisationsdrucktopf to the rim with water, replace and lock the lid, turn the red switch to the closed (ZU) "0" position and than pressurize. The safety valve must open between 2.6-2.8 bar. If the valve does not open, change it or return the unit to the manufacturer. In accordance with safety regulations only the complete unit with matching lid and body should be checked and, if necessary, repaired by the manufacturer with original spare parts. To prolong the life of the unit a regular cleaning outside and inside is recommended. The water should be changed several times daily, depending on frequency of use, since it enriches itself with aggressive components from models and thus possibly can attack the resistant anti-corrosion coating.

We recommend that the periods for retesting registered in the Declaration of Conformity of the performed pressure and inspection test be adhered to in your own interest.

EC-Declaration of Conformity

Description of the unit: DC Polymerisationsdrucktopf (different item no. see above)

Type: Pressure vessel for polymerisation

EC-Guidelines: Pressure Equipment Directive RL 97/23/EC

Capacity: 4.0 l

Maximum pressure: 2.5 bar

Operating pressure: 2.0 – 2.2 bar

Checking pressure: 4.0 bar

Compressor supply: 8.0 bar

Working temperature: max. 50° C / 122° F

Production material: Aluminium / Inside coating Teflon

The serial no. is engraved in handles of lid and body of the unit.

DC Dental Central declares that above described unit corresponds to the relevant fundamental safety and health requirements of the EC-Guidelines due to its conception and design.

This declaration become invalid if the unit is modified or altered in any way without our prior consent. The unit is equipped with a manometer with marking (red line) for the maximum pressure and

a safety valve, complying with TÜV requirements for components, set at 2.5 bar. Recurrent testing of pressure vessel must be ordered by user.

We recommend to test the unit at least every 5 years. The unit must only be operated in accordance with the operating instructions.

Warranty

Item no. 714440 DC Polymerisationsdrucktopf 4 l. with inflating hose

Serial No.: _____

Item no. 714442 DC Polymerisationsdrucktopf 4 l. without inflating hose

Serial No.: _____

Our devices are the result of careful manufacturing and are subject to constant quality controls. If nevertheless a reason should exist for complaint, please consider the following terms of warranty. DC Dental Central grants a warranty for a period of 12 months from the date of the original purchase. It covers defects arising provable from faults in material or manufacturing. This warranty is not transferable.

Warranty performance neither causes an extension of the warranty period nor a further warranty for replacements. Claims of a more far-reaching nature, in particular such on compensation are impossible. We are also not responsible for consequential damages.

ATTENTION: Damages resulting from failure to observe the operating instructions, from negligent handling or misuse, from action by unauthorised persons or by an inevitable event are excluded from the warranty performance.

The performance of the warranty service must be carried out exclusively by the manufacturer and the completion by the responsible dealer.

Costs for shipping and packaging are in charge of the buyer. Transport damages are not subject to warranty claim and must be reported immediately.

Name of client: _____

Date of purchase: _____

Stamp of dealer