**REF 701025**  
transparent /  
transparency

**REF 701027**  
rosa / pink

Inhalt / Content:

**50** Oberkieferplatten  
*Plates maxilla*

**LOT**

25° C  
77° F

Stand 01/2022  
Made in Germany

 DC DentalCentral GmbH • Owidenfeldstr. 6  
D-30559 Hannover • Germany • [www.dental-central.de](http://www.dental-central.de)

**DC DENTAL CENTRAL**

# DC Tray

Lichthärtendes Löffelmaterial  
*Light-curing tray material*

## D Gebrauchsinformation

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## GB Instructions for use

The edge course for the individual tray is marked on the plaster model. Attach space-holders if required. If no space-holders are needed, insulate the plaster model with DC Separating.

Remove protective film and place the DC Tray on the model, shape it and cut at the edges as marked. The transparency of the DC Tray transparent allows very accurate shaping. Excess material can easily be used to make the handle.

DC Tray is polymerised in a normal light device (blue or UVA light). Illuminate the shaped spoon for about 5 minutes from behind.

The shaped tray is ready for use and further processing once it has been polymerised. Any slight inhibition layer produced by oxygen in the air can be removed by wiping the tray with a cloth dipped in DC Clean or an alcohol-based detergent.

### Storage

Store in a dark place

### Precautions

In its unhardened state DC Tray contains methacryl acid ester. Operators who are allergic or hypersensitive should use appropriate protection. It is essential that polymerisation is sufficient and complete, as illumination for an insufficient length of time or at insufficient intensity leads to a loss of strength.