



**DC DENTAL CENTRAL**

# DC Trays

# DC DENTAL CENTRAL

# DC Trays

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## **GB** Disinfection, care and sterilisation of instruments

### **Informations about aluminum-trays**

- **DON'T USE ULTRASOUND** – risks to deteriorate the anodizing and the color
- **DON'T USE THERMODISINFECTOR**
- **ATTENTION:** some products of decontamination are **NOT COMPATIBLE** with aluminium: follow always the enclosed instructions.
- Follow the enclosed instructions of the products you are going to use for the cleaning and remember the incompatibility with aluminium.
- Don't use alkaline products, some products too acid, or products with Sodium or Potassium.
- The PH from 4 to 8 doesn't produce any particular risk.
- The following way of sterilization can be used in safety, always following the enclosed instructions of the producer:

Autoclave: 134° C

Chemiclave: 132° C

Dry Heat: 180° C

## General instructions

*Sterilisation is no substitute for cleaning! If an instrument is exposed to high temperatures before being cleaned first, this can cause stains to be burnt into the instrument.*

Before manual cleaning, the instruments must be disinfected by suitable disinfectants which should not contain any phenol if possible.

Instruments should always be cleaned in a noncorrosive, natural cleaning agent with little foaming.

Stubborn impurities should be removed with a soft brush (never with steel wool, drill brushes or similar). The instrument should then be rinsed in distilled water. Many dental surgeries use ultrasonic units or the thermodisinfectant: here again, care is required to ensure that the cleaning agents and chemicals being used do not damage the metal. And always remember: these methods clean but do not sterilise!

No matter which sterilisation method you use, you should always check your equipment first for any dirt remains and organic or mineral deposits. Both can be transferred to the instruments and cause corrosion. Never mix stainless steel with aluminium, brass or copper.

This could cause electrolysis and also produce etching or corrosion.

Only use distilled or demineralised water. High mineral levels or water that is too soft can cause permanent stains.

Whether you use the autoclave, dry heat or germicidal solvents, you should always follow the manufacturer's instructions precisely when it comes to the temperatures and times involved.

Careful drying is extremely important. Any remaining water or condensation can cause rust or corrosion to an instrument. This applies particularly when packed sterilised, or when the autoclave was opened prematurely.

It is advisable to remove any remaining moisture with a sterile cloth. Instruments consisting of several parts, e.g. pliers, gouging pliers, scissors etc. should always be well oiled. Hinged instruments should always be opened for cleaning and sterilisation.

Although stainless steel has outstanding corrosion resistance, contact with certain chemicals should still be avoided.

The following chemicals should be avoided completely: Aqua Regia, iron chloride, sulphuric acid, hydrochloric acid, iodine.