D Gebrauchsinformation

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GB Instructions for Use

DC Etch+Bond 1 is a light-curing, single-component, self-etching bonding agent for use in the adhesive restorative technique in combination with light-curing composites and compomers.

DC Etch+Bond 1 consists of mono- and bifunctional amide derivatives of methacrylic- and acrylic acid comprising phosphate, carboxylic acid and hydroxyalkyl moieties, water, dispersed silicon dioxide, initiators and stabilizers.

Indication

Establishment of a sound bond on ground enamel and dentin surfaces in direct restorative procedures with light-curing composites and compomers.

Contraindications

- Direct pulp capping
- In combination with dual- or self-curing composites
- If the patient is allergic to any of the ingredients in DC Etch+Bond 1

Interactions

Phenolic and other substances (e.g. zinc-oxide eugenol) may inhibit polymerization. They must not come in contact with DC Etch+Bond 1.

Side effects

In rare cases, repeated contact may lead to a sensitization to monomers used. In these cases, the material should no longer be used.

Application

1. The cavity is prepared according to the principles of the adhesive restorative technique and the treatment field is isolated.
2. DC Etch+Bond 1 is applied with a brush and is massaged into the enamel and dentin surfaces for 30 s. Subsequently, it is dispersed with a strong stream of air in order to ensure evaporation of the water solvent. No mobile liquid film should be visible.
3. The layer is light-cured for 10 s (min. 500 mW/cm²).
4. Subsequently, the restoration is completed as indicated in the directions for use of the composite/compomer employed.

Warnings

- DC Etch+Bond 1 is an irritant. Contact with the skin, mucous membrane and eyes should be avoided. In case of contact, skin/mucous membranes should be rinsed with water. In case of contact with the eyes, they should be immediately rinsed with copious amounts of water and an ophthalmologist should be consulted.
- Commercial medical gloves do not provide protection against the sensitization to methacrylic acid derivatives. Thus contact should also be avoided when wearing gloves.
- Storage at 2 - 8° C / 36 - 46° F
- Expiration date: see imprint on bottle/packaging
- Close bottle immediately after use. Exposure to light will cause premature polymerization.
- Store out of the reach of children!
- For dental use only!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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