DC Alginat Color

D Gebrauchsinformation

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GB Instruction for use

DC Alginat Color is a lead-free alginate showing three different colours during the course of setting. The sequence of colours guarantees accuracy during impression taking even for the unpracticed user. It is germicidal and thus facilitates impression's hygiene. A particular anti-bacterial additive disinfects the impression after removing it from the mouth. Thus, dimensional variations due to immersion into disinfectant solutions are eliminated. Suitable for all impression techniques, e.g. prosthetics, work models or orthodontic models.

Mixing ratio:	17 g powder : 39 ml water
Mixing time manually:	approx. 45-50 sec.
Processing time:	approx. 90 sec.
Setting time:	approx. 120 sec.

(from start of mixing) at a water temperature of 23° C.
Cold water retards, warm water accelerates setting.

Working method:

Invert the can a few times before using in order to loosen the powder.
A measuring set facilitates a precise measuring:

- 1 measuring spoon = 17 g
- 1 measuring cup = 39 ml

Fill the powder into a clean, dry and flexible mixing bowl.
Add water and mix with a spatula until a creamy, pink, soft mixture is obtained.

Impression taking:

We recommend to use a perforated impression tray. Before, the patient should rinse the mouth with warm water. The alginate should be brought into the mouth when the colour turns to white. It is to be held using slight pressure for approx. 30 sec. After having removed the impression tray, the impression is to be rinsed under cold running water.

Important:

Impression should be cast immediately after removal from the mouth.
If this is not possible, the impression is to be covered by a damp cloth and to be stored in a polyethylene bag to keep dimensional variations as small as possible.

Storage:

Closed, room temperature and dry.

ATTENTION:

The impression contains silicogeneous fine dust. The dust formation has been reduced, however breathing dust should be avoided! Gloves should be worn.