DC Alginat Quick

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Str. 2 • D-22946 Trittau • Germany • www.dental-central.de

D Gebrauchsinformation

Vorteile:

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GB Instruction for use

Advantages:

- not injurious to health: no troublesome dust
- aseptic
- clean handling
- rapidly setting
- easy detaching from jaw and teeth
- constant flexible precise impression
- optimum definition of details
- smooth surface of the model

Mixing ratio: 17 g powder: 39 ml water

Mixing time manually: approx. 30 sec.

Processing time: approx. 45-60 sec.

Setting time: approx. 90 sec.

(from start of mixing) at a water temperature of 23° C.

Cold water retards, warm water accelerates setting.

Working method:

Invert the can a few times before using in order to loosen the powder. A measuring set facilitates a precise measuring:

1 measuring spoon = 17 g

1 measuring cup = 39 ml

Strew the powder into the appropriate quantity of water and spatulate thoroughly for a homogenous paste. The perforated impression tray is filled with the mix and inserted into the oral cavity. During the period of intraoral impression taking the impression tray should be supported slightly. After having removed the impression tray, the impression is to be rinsed under cold running water.

Important:

Impression should be cast immediately after removal from the mouth. If this is not possible, the impression is to be covered by a damp cloth and to be stored in a polyethylene bag to keep dimensional variations as small as possible.

Storage:

Closed, room temperature and dry.

ATTENTION:

The impression contains silicogeneous fine dust. The dust formation has been reduced, however breathing in of dust should be avoided!