DC Keramik Ätzgel

Stand 12 / 2010

DC Dental Central Großhandelsges. mbH • Carl-Zeiss-Str. 2 • D-22946 Trittau • Germany • www.dental-central.de • E-Mail: office@dental-central.de

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Allgemeines:

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GB Instructions for use

Areas of application:

All-porcelain restorations, i.e. onlays, inlays, crowns or veneers.

Application:

1. Apply DC ceramic etching gel to the porcelain surface with a small disposable brush or disposable application tip. The recommended reaction time for DC ceramic etching gel is 60 secs.
2. Acid residue should be completely removed after the reaction time. The surface can be steam cleaned, rinsed with plenty of water or cleaned with distilled water in an ultrasonic cleaner to remove the acid residue.
3. Dry thoroughly. Then apply bonding agent etc. according to the manufacturer's instructions.

Important:

The porcelain surface should be etched immediately before fitting or after try-in. Avoid contaminating the etched surfaces of the porcelain (do not brush).

If the surface is contaminated after etching, it should be briefly re-etched (approx. 10 secs.). Etching with DC ceramic etching gel is not required with other bonding systems, e.g. with aluminium oxide porcelain (adhesive technique).

Content: 6 ml

Contains: < 5 % hydrofluoric acid

- Only for use in Dentistry!
- Syringe must be immediately closed after use.
- Keep out of reach of children.
- Store at max. temperature at 25° C.

Safety instructions:

- R26/27/28 Very toxic on inhalation, in contact with skin and if swallowed
- R35 Causes severe burns
- S7/9 Keep container tightly closed and in a well-ventilated place
- S26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice
- S28 After contact with skin wash immediately with plenty of water
- S36/37/39 Wear suitable protective clothing, gloves and eye/face protection
- S45 In case of accident or if you feel unwell seek medical advice immediately (show the label where possible)
- S60 This material and/or its container must be disposed of as hazardous waste

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