DC Keramik Ätzgel

DC DentalCentral GmbH • Owienfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

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Allgemeines:

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GB Instructions for use

Areas fo application:

All-porcelain restorations, i.e. onlays, inlays, crowns or veneers.

Application:

1. Apply DC ceramic etching gel to the porcelain surface with a small disposable brush or disposable application tip. The recommended reaction time for DC ceramic etching gel is 60 secs.
2. Acid residue should be completely removed after the reaction time. The surface can be steam cleaned, rinsed with plenty of water or cleaned with distilled water in an ultrasonic cleaner to remove the acid residue.

3. Dry thoroughly. Then apply bonding agent etc. according to the manufacturer's instructions.

Important:

The porcelain surface should be etched immediately before fitting or after try-in. Avoid contaminating the etched surfaces of the porcelain (do not brush).

If the surface ist contaminated after etching, it should be briefly re-etched (approx. 10 secs.). Etching with DC Ceramic Etching Gel ist not required with other bonding systems, e.g. with aluminium oxide porcelain (adhesive technique).

Safety instructions:

H300 *Fatal if swallowed.*

H310 *Fatal in contact with skin.*

H330 *Fatal if inhaled.*

H314 *Causes severe skin burns and eye damage.*

- Only for use in Dentistry!
- Syringe must be immediately closed after use.
- Keep out of reach of children.
- Store at max. temperature at 25 °C.

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