# DC DENTAL CENTRAL

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# Instructions for use and reprocessing of DC products for endodontics

DC DentalCentral GmbH • Owidenfeldstr. 6 • D-30559 Hannover • Germany • www.dental-central.de • E-Mail: office@dental-central.de

## GB Instructions on reprocessing acc. to DIN EN ISO 17664

### General basics

All instruments must be cleaned, disinfected and sterilized prior to each use; this also applies to first use with instruments delivered non-sterile. Effective cleaning and disinfection is an indispensable prerequisite for effective sterilization. Special instructions on cleaning/sterilization are to be taken from the directions for use. In addition, the operating instructions of your practice devices must be observed. Within the framework of your responsibility for the sterility of the instruments during use, please always observe that only sufficient device-specific and product-specific validated methods are used for cleaning/disinfection and sterilization, that the used devices (disinfector, sterilizer) are regularly maintained and checked and that the validated parameters are complied with in each cycle. **Furthermore, please observe all applicable legal provisions as well as the rules of hygiene of the medical practice.**

This applies especially to the different specifications regarding an effective inactivation of prions. For your own safety, always wear gloves if you are handling with contaminated instruments.

### Cleaning and disinfection

#### Basics

If possible, a mechanical method (disinfector) should be used to clean and disinfect the instruments. Due to its clearly lower efficacy and reproducibility, a manual method – also using an ultrasonic bath – should only be used if no mechanical method is available.

The pretreatment is to be performed in both cases.

#### Pretreatment

Directly after use (within a max. of 2 hours), rough impurities must be removed from the products.

Never let pulp and dentine residues dry up! Directly after use in a patient, put the instruments into the DC Endo Clean Stand (storage period of max. 2 hours) filled with a suitable cleaning agent / disinfectant for intermediate storage and pre-disinfection / cleaning. A clean DC Endo Clean Stand with new foam disk must be used for each patient.

Afterwards, clean the instruments under running water or in a disinfectant solution in order to remove impurities; the disinfectant should be aldehyde-free (otherwise fixation of blood impurities may occur). Only use disinfectants whose efficacy was officially tested (e.g. DGHM\*\* / VAH\*\*\* listing, CE labeling or FDA approval) and which are suitable for the disinfection of instruments and compatible with these instruments.

For a manual removal of impurities, only use a clean soft brush or a clean soft cloth which you only use for this purpose, but never metal brushes or steel wool. Please note that the disinfectant used for pretreatment only serves personal protection and cannot replace the disinfecting step to be performed later on after the cleaning is completed.

### MECHANICAL CLEANING / DISINFECTION

The instruments should not show any visible impurities. Instruments which are clearly damaged (bent, contorted, broken) must be disposed of immediately.

#### – Thermal disinfection (disinfector/RDG\*)

When choosing the disinfector, please make sure that,

- the disinfector was tested for efficacy according to DIN EN ISO 15883,
- if possible, a tested program for thermal disinfection (min. 10 minutes at 93° C or A0 value > 3000) is used,
- the program used is suitable for the instruments and contains sufficient backwash cycles,
- only sterile or germ-free as well as endotoxin-free water (e.g. high purified water HPW) is used for rinsing and that
- the disinfector is regularly maintained and checked.

When choosing the cleaning agent system to be used,

please make sure that

- this system is generally suitable to clean these instruments,
- as long as no thermal disinfection is used, a disinfectant officially tested for efficacy (e.g. DGHM\*\* / VAH\*\*\* listing, CE labeling or FDA approval) is used in addition and is compatible with the used cleaning agent and that
- the used chemicals are compatible with the instruments (see chapter "material resistance").

The concentrations specified by the manufacturer of the cleaning agent and, if

applicable, the disinfectant must imperatively be observed.

#### Procedure:

1. Put the instruments into a kit, a rack or a container in order to prevent any contact between them.
2. Insert the instruments in the disinfector.
3. Start the program.
4. After the program is terminated, remove the instruments from the disinfector.
5. Check and wrap the instruments as immediately as possible after removal. Furthermore, please ensure that there is no contact between instruments / products.

### MANUAL CLEANING AND DISINFECTION

The instruments should not show any visible impurities. Instruments which are clearly damaged (bent, contorted, broken) must be disposed of immediately.

When choosing the cleaning agents and disinfectants to be used, please make sure that

- these are generally suitable to clean or disinfect instruments,
- the cleaning agent, if applicable, is suitable for ultrasonic cleaning (no foaming),
- a disinfectant officially tested for efficacy (e.g. DGHM\*\* / VAH\*\*\* listing, CE labeling or FDA approval) is used and is compatible with the used cleaning agent,
- the used chemicals are compatible with the instruments (see chapter "material resistance").

Combined cleaning agents / disinfectants should only be used in case of very low pre-load of the instruments (no visible impurities). The concentrations and contact times specified by the manufacturer of the cleaning agents and disinfectants must imperatively be observed. Only use freshly prepared solutions, only sterile or germ-free as well as endotoxin-free water (e.g. Aqua purificata (PW)) and only filtered air for drying.

#### Procedure:

##### 1. Cleaning

- Put the instruments into a kit, a rack or a container in order to prevent any contact between them.
- Place the instruments horizontally into the cleaning bath for the specified contact time so that they are sufficiently covered (if appl. supported by ultrasound or cautious brushing with a soft brush).
- Afterwards, remove the instruments from the cleaning bath and rinse them thoroughly for at least 1 minute with water.

##### 2. Disinfection

- Place the cleaned and checked instruments into a disinfectant bath in a disinfection box for the specified contact time with the instruments being sufficiently covered.
- Afterwards, remove the instruments from the disinfectant bath and rinse them thoroughly for at least 1 minute with water.
- Check, dry and wrap the instruments as immediately as possible after removal.

Furthermore, please ensure that there is no contact between instruments.

#### Check

Check all instruments after cleaning or cleaning / disinfection. Instruments are to be sorted out immediately if they show deficiencies, such as:

- plastically deformed
- instrument is bent
- windings are loose
- cutting surfaces are damaged
- edges are blunt
- strength labeling is missing
- corrosion

For information on numeric restriction of reuse please refer to "reusability". Instruments which are still soiled must be cleaned and disinfected again.

#### Wrapping

Please wrap the instruments in the Endo sterilization trays and then in disposable sterilization wrappings (single wrapping) which meet the following requirements:

- according to DIN EN 868/ANSI AAMI ISO 11607
- suitable for steam sterilization (temperature resistance up to at least 141° C (286° F) sufficient vapor permeability)

### Sterilization

Only the subsequently mentioned sterilization methods are to be used for the sterilization; other sterilization methods are not allowed.

#### Steam sterilization

- fractionated vacuum method or gravitation method<sup>1</sup> (with sufficient product drying)
- steam sterilizer according to DIN EN 13060 or DIN EN 285
- validated according to DIN EN 554/ANSI AAMI ISO 11134 (valid picking and product-specific performance evaluation)
- maximum sterilization temperature 138° C (280° F); plus tolerance according to DIN EN 554/ ANSI AAMI ISO 11134)
- sterilization time (exposure time at sterilization temperature) min. 20 minutes (at 121° C (250° F) or 5 minutes<sup>2</sup> at 132° C / 134° C (270° F)

The lightning sterilization method or the sterilization of unwrapped instruments is by no means allowed.

Moreover, do not use any hot-air sterilization, no radiation sterilization, no formaldehyde or ethylene oxide sterilization, as well as no plasma sterilization.

#### Storage

After sterilization, the instruments must be stored in the sterilization packaging in a dry and dust-free place.

#### Material resistance

When choosing the cleaning agents and disinfectants please make sure that no phenolic

Product	Material	Special indications on cleaning/ sterilization	Reusability	Possible damage
DC Hedstroem files	Stainless steel and temperature-resistant plastic		Cleaned and undamaged instruments can be reused 8 to 10 times	Ruptures at the plastic grip, corrosion at the working part
DC NiTi K files DC NiTi reamer DC K Reamer	NiTi alloy and temperature-resistant plastic, stainless steel		Cleaned and undamaged instruments can be reused up to 7 times	Ruptures at the plastic grip, corrosion at the working part
DC NiTi spreader	NiTi alloy and temperature-resistant plastic		Cleaned and undamaged files can be reused without limit; please watch out for any wear and tear	Ruptures at the plastic grip, corrosion at the working part
DC paste carrier	Stainless steel		Cleaned and undamaged instruments can be reused up to 5 times	Corrosion at the working part shaft
DC Endo storage, DC Endo meter, DC NiTi clean stand DC Endo file measure set	Plastic or aluminum	For sterilization, the product must be welded in a disposable sterilization wrapping		
DC Foam inserts for DC Endo clean stand + Twin Block	Foam	Can be autoclaved prior to single use	For single use only	
Silicone stopper	Silicone	Silicone stopper should be removed prior to cleaning/disinfection and be separately cleaned/disinfected	Single use of the stoppers is recommended	
DC gutta-percha points color	Gutta-percha, zinc oxide, barium sulfate	Cold disinfection, e.g. in med. alcohol	For single use only	Deformation
DC paper points color		Cold disinfection, e.g. in med. alcohol	For single use only	Deformation

<sup>1</sup> Using the less effective gravitation method is only allowed if the fractionated vacuum method is not available.

<sup>2</sup> or 18 minutes (inactivation of prions)

## GB Instruction for use

We recommend to always work with a rubber dam!

### 1. Hand instruments

Product	Application	Specific details
DC Hedstroem-Files	Pushing and pulling motion without rotation. Rotating motion would cause the sharp blades to block. If the canal was enlarged with Reamers or K-files the following Hedstroem file should be one size smaller or equal.	
DC NiTi K-Files	<b>NiTi K-Files must always be used manually!</b> Pushing an pulling motion without rotation. Rotating motion would cause the sharp blades to block.	NiTi-Instruments contain nickel and should not be used for individuals with known allergic sensitivity to this metal.
DC NiTi Spreader	Lateral condensation of gutta-percha points. The spreader is placed between the points and pushed carefully in apical direction.	NiTi-Instruments contain nickel and should not be used for individuals with known allergic sensitivity to this metal.
DC NiTi Reamer DC K Reamer	Pushing/rotating motion (reaming motion) max. 90° clockwise.	NiTi-Instruments contain nickel and should not be used for individuals with known allergic sensitivity to this metal.

### 2. Instrument for handpiece

Product	Application	Specific details
DC paste carrier	With the handpiece shut off, submerge paste carrier in filler material, insert paste carrier close to the apex, and with max. 400-800 rpm rotate the filler material into the canal whilst slowly withdrawing the instrument from the canal.	

### 2. Endo-boxes

Product	Application
DC Endo-boxes	Sort the instruments into the module. Close the wash box. Both parts must lock into place. Then either shrink wrap and autoclave or remove the module with instruments and autoclave them in an endo-box (with autoclave paper inserted).
DC Endo Clean stand	Fill the container with a suitable disinfection solution (e.g. NaOCl), cover with DC foam inserts and fasten with the collar. During treatment stick the instruments in the foam inserts for interim storage and cleaning.

\* cleaning and disinfectant devices

\*\* German Association for Hygiene and Microbiology,

\*\*\* German Association for Applied Hygiene