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GB User information and reprocessing

All instruments are developed and designed for a specific application. Improper use can lead to tissue injuries, premature wear, destruction of the instruments and danger to the user, the patient or a third party.

Rotary instruments for medical applications may only be employed by doctors or by other individuals who are familiar with the use of such instruments as a result of their education or training and their particular experience.

Proper application

Care must be taken that only technically and hygienically serviced and cleaned turbines, hand pieces and contraangles are employed.

- Insert instruments as deeply as possible.
- Before placing the instruments on the object, be sure that the desired speed has been reached.
- Avoid tilting or leveraging because of an increased danger of breakage.
- Depending on the particular application, the use of protective goggles may be recommended.
- Users must avoid handling the instruments without wearing protective gloves.
- Improper application leads to poor working results and increases the risk of injury.

Recommended speed

- Failure to observe the recommended speed may lead to increased safety risk.
- When exceeding the maximum recommended speed, long and pointed instruments tend to vibrate, which may lead to instrument damage.
- Strong centrifugal forces can appear when employing tools with a working diameter larger than the diameter of the shaft, which may lead to bending of the shaft and/or to the breaking of the instrument. The maximum recommended speed should by no means be exceeded.

Contact pressure

Avoid excessive contact pressure since this will lead to defective blades when employing cutting instruments. At the same time, this may also lead to increased heat.

- Excessive contact pressure in grinding instruments leads to excessive wear and/or to an excessive generation of heat.
 - Through overheating, excessive contact pressure can lead to damage of the dental pulp.
- Fractured blades result in undesirably rough surfaces. The breakage of an instrument cannot be excluded.

Cooling

- In order to avoid the undesirable development of heat, it is mandatory to cool sufficiently using an air/water spray (with a minimum of 50 ml/min).
- For FG instruments, with a total length of over 22 mm or a head diameter of over 2 mm, additional external cooling is necessary.
- Insufficient water-cooling leads to irreversible damage of the tooth and the surrounding tissues.

Discarding of tools

- Fractured and irregular blades lead to vibrations. High contact pressure leads to fracturing of the preparation edges and to rough surfaces.
- Bent or worn instruments, or instruments which do not run smoothly, must be discarded immediately.

Cleaning and hygiene

For disinfection and sterilization, care must be taken to see that the procedure being carried out is suitable for the respective instrument.

- Before using a tool on a patient for the first time, and immediately after each use, rotaries must be disinfected, cleaned and sterilized. Until first used, they may be stored in the original package at room temperature, free from dust and moisture.
- The storage of rotaries should occur in hygienically-serviced stands, bowls or in other suitable containers.
- All packaged instruments have not been sterilized and must consequently be sterilized before use on the patient. The sterilization is to be carried out using a validated process. Cleaning and sterilization of the instruments which is not suitable/not carried out according to instructions may lead to an infection in the patient involving dangerous pathogens.
- For instruments which are not protected from corrosion, disinfectants and cleaning agents with corrosion prevention must be employed.
- All suitable materials may be sterilized in an autoclave with damp heat according to EN DIN 554 (e.g., at 134° C for a period of 3 min).
- Contact with H₂O₂ (hydrogen peroxide) must be avoided since this will cause damage to carbide components.
- Avoid temperatures over 180° C. Exceeding this temperature leads to a reduction in the lifetime of the instrument.
- Carbide rotary instruments and instruments which are not corrosion-resistant corrode in the thermodesinfector. This leads to discolorations and a reduced lifetime.
- Any handling of contaminated instruments must be carried out with suitable protective gloves.
- All parts must be cleaned and disinfected sufficiently prior to sterilization. In the event of manual cleaning, blood and tissue parts must be removed from the instruments using a suitable brush.
- Directions for use, duration of exposure and suitability of disinfectants and cleaning substances for certain types of instruments are to be taken from the instructions of the manufacturers of these products.

Storage

- Store instruments in a dry place.
- Protect instruments from sunlight, heat and dust.