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GB User information and reprocessing

By purchasing these instruments you obtain a high-quality product whose proper handling and usage are described below. To minimize danger for patients and users, we ask you to carefully comply with these instructions for use. Use, disinfection, cleaning and sterilization of the instruments may be performed by instructed specialists only.

Control

Prior to any use, the instruments must be examined for integrity. Surface damage such as scratches, flaws, jags, dents etc. as well as deformed parts mean that the instruments must not be used. In this case, the products are to be fed into the normal disposal process of the hospital. Do not use any damaged products!

Field of application

Stomatoscope for dental examinations.

Handling

The instruments must not be used for any other than the intended purpose since this may result in damage to or breaking of parts of the instruments.

Risks

None known.

Complications

Generally there are no complications.

Safety precautions

It is pointed out that in case of convulsions or seizures (e. g. epilepsy) the mirror in the oral cavity may cause injury (fragment injury). As a safety precaution, mouth gags are to be used if required.

Disposal

If it should not be possible any more to reprocess the instruments, they are to be fed into the normal disposal process of the hospital.

Materials

The materials used are stainless steels according to 1.4016 and glass.

PROCESSING INSTRUCTIONS

Procedure:

- Cleaning
- Disinfection
- Pressurized vapour sterilization in the autoclave (DIN EN ISO 17665-1)

Attention: Dry heat sterilization no longer conforms to the current state of the art.

Warning:

The instruments are shipped non-sterile and must be cleaned and, if necessary, disinfected and sterilized before use.

The instruments may be processed only by persons with the necessary training and skills and are capable of judging the potential risks and their respective effects.

Limitations on reprocessing:

Frequent proper reprocessing has little impact on the instruments. Normally, the end of product lifetime is determined by

wear and damage due to use. In this case, the instruments are to be fed into the normal disposal process of the hospital. Do not use any damaged products!

DC endo-mouth mirrors: Handle the reflecting surface with care, because a razor-thin coating is coated on it.

INSTRUCTIONS:

Place of use:

Remove surface contaminations using a single use towel/ paper towel. It is recommended to reprocess the instruments as early as possible after use.

Immediately after use, manual disinfection may be performed to reduce infection risk for the users. To this end, the instruments are immersed into a disinfectant with demonstrated efficacy (e. g. Certificate by DGHM/VAH). Make sure that the instruments are immersed into the disinfectant solution completely and without bubble formation. Pay attention to the instructions by the disinfectant manufacturer relating to concentration, residence time and temperature.

Preparation for decontamination:

If the instruments can be disassembled, disassemble them prior to reprocessing.

CLEANING: MANUALLY

- Immediately after use, clean the outside of the instruments using a soft towel, sponge and suitable soft brushes in a cleaning solution.
- The cleaning solution used must be suitable for cleaning steel products and glass.
- The instruments must be completely covered in the solution. Do not overload the container.
- Precisely comply with concentration, residence time and temperature according to the manufacturer's instructions.
- The cleaning solution should be a non-foaming solution suitable for manual cleaning.
- Carefully rinse the interior of instruments with channels or cavities.
- Perform brushing and all further cleaning steps below the surface to avoid squirting of contaminated liquid.
- Take the instruments out of the cleaning solution.
- Carefully dry after rinsing with clear, running water.
- Change the cleaning solution at least once daily or immediately after visible contamination
- For persistent contamination, e. g. with phosphate cements, we recommend «DC Orange» by Dental Central.

Ultrasound cleaning

Ultrasound cleaning is not recommended.

DISINFECTION: MANUALLY

- The disinfectant with demonstrated efficacy (e. g. Certificate by DGHM/VAH) to be used must be suitable for disinfection of glass and steel products (e. g. «DC Instrumentenbad»).
- Immerse cleaned instruments into a disinfectant trough.
- Fill all channels and cavities with disinfectant solution so that there are no bubbles formed.
- Cover the trough with appropriate lid.
- Precisely comply with concentration, residence time and temperature according to the manufacturer's instructions. Wrong concentration and too long residence times result in damage to the material.
- Wearing fresh single-use gloves, take the instruments out of the disinfectant solution.

Neutralisation/ rinsing

- Immerse disinfected instruments into bin/trough with micro-biologically acceptable/ sterile water; use fresh water for each instrument.
- Thoroughly rinse the outside of the instruments and all channels and cavities to remove residual disinfectant.
- Take the instruments out of the water.

Drying and functional test

- Dry external surfaces with lint-free cloth and pressurized air.
- Completely dry all channels and cavities using pressurized air.
- Test instrument for proper function.

Maintenance, control

Separate out damaged instruments and dispose of them. Prior to disposal, it is mandatory to sterilize the instruments.

Packaging:

Sterilization and storage in foil is not possible with stomatoscopes.

STERILIZATION:

Recommended method for sterilization:

Vapour sterilization with saturated vapour and fractionated vacuum (DIN EN ISO 17665-1)

Recommended temperature:

134° C

Recommended pressure:

3 bar

Holding time:

≥ 5 min

Drying time:

≥ 15 min

Alternative method:

Dry heat sterilization at # maximally 180° C.

Attention: Dry heat sterilization no longer conforms to the current state of the art.

Storage:

Store sterilized instruments in sterile material packaging in a closed locker, protected from dust, moisture and temperature variations.

Additional information:

Further information on reprocessing of medical products:

- Internet: <http://www.rki.de>
- Internet: <http://www.a-k-i.org> brochure «Instrumentenaufbereitung in der Zahnarztpraxis richtig gemacht»
- Hygiene requirements in reprocessing medical products: Recommendations by the Committee for Hospital Hygiene and Infection Prevention at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Products (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM)
- Information on reprocessing of sterilisable medical products DIN EN ISO 17664

The instructions listed above were validated by the manufacturer of medical products as SUITABLE for preparing a medical product for reuse. The person performing the reprocessing shall be responsible for obtaining the desired results of the actually performed reprocessing with the equipment, materials and persons in the reprocessing system. To this end, validation and routine monitoring of the process will normally be required.

Warranty:

The products are manufactured from high-quality materials and undergo quality control prior to shipping. If there should nevertheless be any flaws, please contact us.

However, we cannot guarantee that the products are suitable for any individual operation. The user shall be responsible for determining this. We shall be exempt from liability for any random or consequential damage.

DC DentalCentral shall be exempt from any liability if these instructions for use have been demonstrably violated.

Attention:

We decline any responsibility for reuse in case any instrument has been used in patients with Creutzfeld-Jakob disease or an HIV infection.